


## Product datasheet for **GA208935**

### Mouse Vwce activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | Vwce CRISPRa kit - CRISPR gene activation of mouse von Willebrand factor C and EGF domains                                                                                                                                                                                                              |
| Symbol:       | Vwce                                                                                                                                                                                                                                                                                                    |
| Components:   | GA208935G1, Vwce gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208935G2, Vwce gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208935G3, Vwce gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q3U515</a>                                                                                                                                                                                                                                                                                  |
| Synonyms:     | 1300015B04Rik                                                                                                                                                                                                                                                                                           |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)