

Product datasheet for **GA208951**

Mouse *Pitpnc1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pitpnc1 CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol transfer protein, cytoplasmic 1
Symbol:	Pitpnc1
UniProt ID:	Q8K4R4
Synonyms:	1110020B03Rik; 5830436L09Rik; A1662802; A1851387; C330017I21Rik; Dnr411; rdgB-beta; RDGGB; RDGGB1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)