

Product datasheet for **GA208982**

Mouse Tex44 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tex44 CRISPRa kit - CRISPR gene activation of mouse testis expressed 44
Symbol:	Tex44
Components:	GA208982G1, Tex44 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208982G2, Tex44 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208982G3, Tex44 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9DA60</u>
Synonyms:	1700019O17Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)