

Product datasheet for **GA209006**

Mouse Tmem106b activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Tmem106b CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 106B |
| Symbol: | Tmem106b |
| UniProt ID: | <u>Q80X71</u> |
| Synonyms: | 2310036D22Rik; 5830455K21Rik; 6430519M21Rik; AI428776; AI661344 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)