

Product datasheet for **GA209010**

Mouse Paqr7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Paqr7 CRISPRa kit - CRISPR gene activation of mouse progesterin and adipoQ receptor family member VII
Symbol:	Paqr7
Components:	GA209010G1, Paqr7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209010G2, Paqr7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209010G3, Paqr7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80ZE4</u>
Synonyms:	2310021M12Rik; AV001994; BB007237; mPR; Mpra; PGLP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)