

Product datasheet for **GA209037**

Mouse Wdr73 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr73 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 73
Symbol:	Wdr73
UniProt ID:	<u>Q9CWR1</u>
Synonyms:	1200011I23Rik; 2410008B13Rik; AI848574; C85352
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)