

Product datasheet for **GA209043**

Mouse Tdrd12 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tdrd12 CRISPRa kit - CRISPR gene activation of mouse tudor domain containing 12
Symbol:	Tdrd12
Components:	GA209043G1, Tdrd12 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209043G2, Tdrd12 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209043G3, Tdrd12 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CWU0</u>
Synonyms:	2410004F06Rik; 2410070K17Rik; ECAT8; EG434165; G1-476-14; repro23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)