

Product datasheet for **GA209048**

Mouse Ddx28 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ddx28 CRISPRa kit - CRISPR gene activation of mouse DEAD (Asp-Glu-Ala-Asp) box polypeptide 28
Symbol:	Ddx28
Components:	GA209048G1, Ddx28 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209048G2, Ddx28 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209048G3, Ddx28 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CWT6</u>
Synonyms:	2410004K13Rik; AI449652; Mddx28
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)