


Product datasheet for **GA209075**

Mouse Ddx42 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ddx42 CRISPRa kit - CRISPR gene activation of mouse DEAD (Asp-Glu-Ala-Asp) box polypeptide 42
Symbol:	Ddx42
Components:	GA209075G1, Ddx42 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209075G2, Ddx42 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209075G3, Ddx42 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q810A7</u>
Synonyms:	1810047H21Rik; AW319508; AW556242; B430002H05Rik; RHELP; RNAHP; SF3b125
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)