

Product datasheet for **GA209108**

Mouse Zbed3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zbed3 CRISPRa kit - CRISPR gene activation of mouse zinc finger, BED type containing 3
Symbol:	Zbed3
Components:	GA209108G1, Zbed3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209108G2, Zbed3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209108G3, Zbed3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D0L1</u>
Synonyms:	0610037K01Rik; 2610005H11Rik; AU018975; AU024588
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)