

Product datasheet for **GA209114**

Mouse *Seh1l* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Seh1l CRISPRa kit - CRISPR gene activation of mouse SEH1-like (<i>S. cerevisiae</i>)
Symbol:	Seh1l
Components:	GA209114G1, Seh1l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209114G2, Seh1l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209114G3, Seh1l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8R2U0
Synonyms:	2610007A16Rik; AW540070; SEC13L; Seh1; SEH1A; SEH1B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)