

## Product datasheet for **GA209138**

### Mouse Shq1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | Shq1 CRISPRa kit - CRISPR gene activation of mouse SHQ1 homolog ( <i>S. cerevisiae</i> )                                                                                                                                                                                                                |
| Symbol:       | Shq1                                                                                                                                                                                                                                                                                                    |
| Components:   | GA209138G1, Shq1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209138G2, Shq1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209138G3, Shq1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q7TMX5</a></u>                                                                                                                                                                                                                                                                           |
| Synonyms:     | 2810403P18Rik; Grim-1                                                                                                                                                                                                                                                                                   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)