

Product datasheet for **GA209148**

Mouse Supt7l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Supt7l CRISPRa kit - CRISPR gene activation of mouse suppressor of Ty 7-like
Symbol:	Supt7l
Components:	GA209148G1, Supt7l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209148G2, Supt7l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209148G3, Supt7l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CZV5</u>
Synonyms:	2610524B01Rik; 6030455L14Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)