

## Product datasheet for **GA209212**

### Mouse Dusp3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Dusp3 CRISPRa kit - CRISPR gene activation of mouse dual specificity phosphatase 3 (vaccinia virus phosphatase VH1-related)
Symbol:	Dusp3
UniProt ID:	<a href="#">Q9D7X3</a>
Synonyms:	2210015003Rik; 5031436003Rik; T-DSP11; VHR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)