


Product datasheet for **GA209223**

Mouse Cdkn3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cdkn3 CRISPRa kit - CRISPR gene activation of mouse cyclin-dependent kinase inhibitor 3
Symbol:	Cdkn3
UniProt ID:	Q810P3
Synonyms:	2410006H10Rik; KAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)