


Product datasheet for **GA209231**

Mouse Wdr44 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr44 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 44
Symbol:	Wdr44
Components:	GA209231G1, Wdr44 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209231G2, Wdr44 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209231G3, Wdr44 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6NVE8</u>
Synonyms:	2610034K17Rik; RAB11BP; RPH11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)