

Product datasheet for **GA209235**

Mouse Dnajc25 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dnajc25 CRISPRa kit - CRISPR gene activation of mouse Dnaj heat shock protein family (Hsp40) member C25
Symbol:	Dnajc25
Components:	GA209235G1, Dnajc25 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209235G2, Dnajc25 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209235G3, Dnajc25 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2ALW5
Synonyms:	2010109C08Rik; 2010203O07Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)