

## Product datasheet for **GA209248**

### Mouse Tmem87b activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Tmem87b CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 87B                                                                                                                                                                                                                                  |
| Symbol:       | Tmem87b                                                                                                                                                                                                                                                                                                          |
| Components:   | GA209248G1, Tmem87b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209248G2, Tmem87b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209248G3, Tmem87b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8BKU8</a></u>                                                                                                                                                                                                                                                                                    |
| Synonyms:     | 2610301K12Rik; 2810431I02Rik; AU014804                                                                                                                                                                                                                                                                           |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)