

Product datasheet for **GA209250**

Mouse Tspyl4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tspyl4 CRISPRa kit - CRISPR gene activation of mouse TSPY-like 4
Symbol:	Tspyl4
Components:	GA209250G1, Tspyl4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209250G2, Tspyl4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209250G3, Tspyl4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VD63</u>
Synonyms:	2610102M01Rik; B230210I21Rik; D10Bwg0791e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)