

Product datasheet for GA209255

Mouse 2610507B1IRik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	2610507B1IRik CRISPRa kit – CRISPR gene activation of mouse RIKEN cDNA 2610507B11 gene
Symbol:	2610507B1IRik
Components:	GA209255G1, 2610507B1IRik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209255G2, 2610507B1IRik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209255G3, 2610507B1IRik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5SYL3</u>
Synonyms:	AU022687; D11Bhm178e; D11Bhm179e; EI; Kiaa0100; mKIAA0100
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

