

Product datasheet for **GA209256**

Mouse **Taf4b** activation kit by **CRISPRa**

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Taf4b CRISPRa kit - CRISPR gene activation of mouse TATA-box binding protein associated factor 4b |
| Symbol: | Taf4b |
| UniProt ID: | G5E8Z2 |
| Synonyms: | 105kDa; 2610524B04Rik; 4932409F03Rik; AW987595; Taf2c2; TAFII-105; TAFII105 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)