

Product datasheet for GA209262

Mouse Wdr43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr43 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 43
Symbol:	Wdr43
Components:	GA209262G1, Wdr43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209262G2, Wdr43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209262G3, Wdr43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6ZQL4
Synonyms:	2610318G08Rik; AU020887; AV024208; mKIAA0007; Wrd43
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)