

Product datasheet for **GA209285**

Mouse Pan3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pan3 CRISPRa kit - CRISPR gene activation of mouse PAN3 poly(A) specific ribonuclease subunit
Symbol:	Pan3
Components:	GA209285G1, Pan3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209285G2, Pan3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209285G3, Pan3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q640Q5</u>
Synonyms:	2700050F09Rik; A430027N15Rik; AU014670
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)