


## Product datasheet for **GA209296**

### Mouse Lins1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Lins1 CRISPRa kit - CRISPR gene activation of mouse lines homolog 1                                                                                                                                                                                                                                        |
| Symbol:       | Lins1                                                                                                                                                                                                                                                                                                      |
| Components:   | GA209296G1, Lins1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209296G2, Lins1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209296G3, Lins1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q3U1D0</a></u>                                                                                                                                                                                                                                                                              |
| Synonyms:     | 2700083B01Rik; AW547219; Lins; Lins2; Wins2                                                                                                                                                                                                                                                                |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)