

Product datasheet for **GA209316**

Mouse HnrnpII activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HnrnpII CRISPRa kit - CRISPR gene activation of mouse heterogeneous nuclear ribonucleoprotein L-like
Symbol:	HnrnpII
Components:	GA209316G1, HnrnpII gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209316G2, HnrnpII gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209316G3, HnrnpII gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q921F4</u>
Synonyms:	2510028H02Rik; 2810036L13Rik; AI256697; AI852082; HnrpII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)