


Product datasheet for **GA209330**

Mouse Tmem161b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem161b CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 161B
Symbol:	Tmem161b
Components:	GA209330G1, Tmem161b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209330G2, Tmem161b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209330G3, Tmem161b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8C2L6
Synonyms:	2810446P07Rik; AI843389
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)