

Product datasheet for **GA209335**

Mouse Tmem135 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem135 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 135
Symbol:	Tmem135
Components:	GA209335G1, Tmem135 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209335G2, Tmem135 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209335G3, Tmem135 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CYV5</u>
Synonyms:	2810439K08Rik; AW319712; PMP52
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)