

Product datasheet for **GA209366**

Mouse Hepacam activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hepacam CRISPRa kit - CRISPR gene activation of mouse hepatocyte cell adhesion molecule
Symbol:	Hepacam
UniProt ID:	<u>Q640R3</u>
Synonyms:	2900042E01Rik; Glialcam; Hepn1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)