

Product datasheet for **GA209368**

Mouse Swi5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Swi5 CRISPRa kit - CRISPR gene activation of mouse SWI5 recombination repair homolog (yeast)
Symbol:	Swi5
Components:	GA209368G1, Swi5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209368G2, Swi5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209368G3, Swi5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8K3D3
Synonyms:	1500019F05Rik; 2900010J23Rik; AI046671
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)