

Product datasheet for **GA209369**

Mouse Ddx41 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ddx41 CRISPRa kit - CRISPR gene activation of mouse DEAD (Asp-Glu-Ala-Asp) box polypeptide 41
Symbol:	Ddx41
UniProt ID:	<u>Q91VN6</u>
Synonyms:	2900024F02Rik; AA958953; ABS; AI324246
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)