


Product datasheet for **GA209380**

Mouse Cnih3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cnih3 CRISPRa kit - CRISPR gene activation of mouse cornichon family AMPA receptor auxiliary protein 3
Symbol:	Cnih3
Components:	GA209380G1, Cnih3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209380G2, Cnih3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209380G3, Cnih3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6ZWS4
Synonyms:	2900075G08Rik; CNIH-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)