


Product datasheet for **GA209408**

Mouse Tmed5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmed5 CRISPRa kit - CRISPR gene activation of mouse transmembrane p24 trafficking protein 5
Symbol:	Tmed5
Components:	GA209408G1, Tmed5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209408G2, Tmed5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209408G3, Tmed5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CXE7</u>
Synonyms:	3110020O18Rik; 4432412D15Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)