


## Product datasheet for **GA209412**

### Mouse Clec4a3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Clec4a3 CRISPRa kit - CRISPR gene activation of mouse C-type lectin domain family 4, member a3
Symbol:	Clec4a3
Components:	GA209412G1, Clec4a3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209412G2, Clec4a3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209412G3, Clec4a3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	3110037K17Rik; Dcir3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)