

Product datasheet for **GA209425**

Mouse Xpot activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Xpot CRISPRa kit - CRISPR gene activation of mouse exportin, tRNA (nuclear export receptor for tRNAs)
Symbol:	Xpot
Components:	GA209425G1, Xpot gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209425G2, Xpot gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209425G3, Xpot gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	1110004L07Rik; 3110065H13Rik; AI452076; C79645; EXPORTIN-T
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)