

Product datasheet for **GA209434**

Mouse Snorc activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Snorc CRISPRa kit - CRISPR gene activation of mouse secondary ossification center associated regulator of chondrocyte maturation
Symbol:	Snorc
Components:	GA209434G1, Snorc gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209434G2, Snorc gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209434G3, Snorc gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CXL7</u>
Synonyms:	Snorc
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)