

Product datasheet for **GA209447**

Mouse Ccdc103 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ccdc103 CRISPRa kit - CRISPR gene activation of mouse coiled-coil domain containing 103
Symbol:	Ccdc103
Components:	GA209447G1, Ccdc103 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209447G2, Ccdc103 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209447G3, Ccdc103 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D9P2</u>
Synonyms:	1700039D13Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)