

Product datasheet for **GA209449**

Mouse 1700031F05Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	1700031F05Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 1700031F05 gene
Symbol:	1700031F05Rik
Components:	GA209449G1, 1700031F05Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209449G2, 1700031F05Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209449G3, 1700031F05Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	RP23-76A20.6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)