


Product datasheet for **GA209464**

Mouse Prss44 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Prss44 CRISPRa kit - CRISPR gene activation of mouse protease, serine 44
Symbol:	Prss44
Components:	GA209464G1, Prss44 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209464G2, Prss44 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209464G3, Prss44 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q402U7</u>
Synonyms:	79H19A; 1700036D21Rik; AU043574; AU043597; TESSP-4; Tessp4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)