

Product datasheet for **GA209468**

Mouse 1700034J05Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	1700034J05Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 1700034J05 gene
Symbol:	1700034J05Rik
Components:	GA209468G1, 1700034J05Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209468G2, 1700034J05Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209468G3, 1700034J05Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80W69</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)