


## Product datasheet for **GA209473**

### Mouse 1700056E22Rik activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                                    |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                                       |
| Description:  | 1700056E22Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 1700056E22 gene                                                                                                                                                                                                                                             |
| Symbol:       | 1700056E22Rik                                                                                                                                                                                                                                                                                                                      |
| Components:   | GA209473G1, 1700056E22Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209473G2, 1700056E22Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209473G3, 1700056E22Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                                            |

#### Product images:



[View online »](#)