

Product datasheet for **GA209476**

Mouse **Tex33** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Tex33 CRISPRa kit - CRISPR gene activation of mouse testis expressed 33
Symbol:	Tex33
Components:	GA209476G1, Tex33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209476G2, Tex33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209476G3, Tex33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D9J2
Synonyms:	1700061J05Rik; Ean57
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)