

Product datasheet for **GA209479**

Mouse Fam177a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Fam177a CRISPRa kit - CRISPR gene activation of mouse family with sequence similarity 177, member A
Symbol:	Fam177a
Components:	GA209479G1, Fam177a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209479G2, Fam177a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209479G3, Fam177a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BR63</u>
Synonyms:	1700047I17Rik; 1700047I17Rik1; C14orf24; D730001C10Rik; Fam177a1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)