


Product datasheet for **GA209480**

Mouse Bpifa3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Bpifa3 CRISPRa kit - CRISPR gene activation of mouse BPI fold containing family A, member 3
Symbol:	Bpifa3
Components:	GA209480G1, Bpifa3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209480G2, Bpifa3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209480G3, Bpifa3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D9J8
Synonyms:	1700058C13Rik; Splunc3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)