

Product datasheet for **GA209485**

Mouse Nipsnap3a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Nipsnap3a CRISPRa kit - CRISPR gene activation of mouse nipsnap homolog 3A
Symbol:	Nipsnap3a
Components:	GA209485G1, Nipsnap3a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209485G2, Nipsnap3a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209485G3, Nipsnap3a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	1700054F22Rik; AI429120; RP23-95M14.6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)