

Product datasheet for **GA209487**

Mouse Tepp activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tepp CRISPRa kit - CRISPR gene activation of mouse testis, prostate and placenta expressed
Symbol:	Tepp
Components:	GA209487G1, Tepp gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209487G2, Tepp gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209487G3, Tepp gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6IMH0</u>
Synonyms:	1700055M20Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)