

## Product datasheet for GA209498

### Mouse Tex35 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tex35 CRISPRa kit - CRISPR gene activation of mouse testis expressed 35
Symbol:	Tex35
UniProt ID:	<u>Q14BK3</u>
Synonyms:	1700057K13Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)