


Product datasheet for **GA209508**

Mouse Ms4a13 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ms4a13 CRISPRa kit - CRISPR gene activation of mouse membrane-spanning 4-domains, subfamily A, member 13
Symbol:	Ms4a13
Components:	GA209508G1, Ms4a13 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209508G2, Ms4a13 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209508G3, Ms4a13 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5FWC3
Synonyms:	1700060E18Rik; Tetm4.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)