

Product datasheet for **GA209527**

Mouse **Pebp4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Pebp4 CRISPRa kit - CRISPR gene activation of mouse phosphatidylethanolamine binding protein 4
Symbol:	Pebp4
Components:	GA209527G1, Pebp4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209527G2, Pebp4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209527G3, Pebp4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>G3X8V3</u>
Synonyms:	PEBP-4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)