

Product datasheet for **GA209575**

Mouse Wdr75 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr75 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 75
Symbol:	Wdr75
Components:	GA209575G1, Wdr75 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209575G2, Wdr75 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209575G3, Wdr75 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3U821</u>
Synonyms:	1300003A18Rik; 2410118I19Rik; C77608
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)