

Product datasheet for **GA209592**

Mouse *Sh3bgrl3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sh3bgrl3 CRISPRa kit - CRISPR gene activation of mouse SH3 domain binding glutamic acid-rich protein-like 3
Symbol:	Sh3bgrl3
Components:	GA209592G1, Sh3bgrl3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209592G2, Sh3bgrl3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209592G3, Sh3bgrl3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q91VW3
Synonyms:	1110004L05Rik; AU022188
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)